



QUICK REFERENCE

OPERATIONS GUIDE


SaltMED™ Inc

206 Birmingham Drive, Cardiff by the Sea, CA 92007
(760) 334.5267 • info@thesaltfacial.com • www.saltmed.com

TABLE OF CONTENTS

SEA SALT RESURFACING	3
AESTHETIC ULTRASOUND	4
LED PHOTOTHERAPY	5 -7
CARE OF THE RECHARGEABLE DEHUMIDIFIER	8
CARE OF THE LED ARTICULATING ARM & PANELS	9
DAILY MAINTENANCE	10
WEEKLY MAINTENANCE	11
SERVICE & SUPPORT	12

HOW TO USE QR CODE



Welcome! This guide utilizes QR code technology. On some cards, you will see a smartphone-camera scannable code that provides an online video or link for more information on the subject of the card. On some older phones you may need to install a QR code app.

POINT IT

Point your phone at the QR CODE

SCAN IT

Take a picture of the QR CODE

ENJOY IT!

Once the QR CODE is received by your phone just relax and enjoy the results

STEP 1: SEA SALT RESURFACING

1. TURN SEA SALT MODULE ON

Press the black On/Off toggle switch on left side of Sea Salt Resurfacing Control Panel to 'ON' (I).


2. SELECT INTENSITY LEVEL

To adjust the Intensity Level: press the round +/- buttons on each side of the Intensity window. Intensity Levels range from 3 (lowest) to 10 (highest).

The system will automatically default to the last Intensity Level used.

3. TURN SALT FLOW ON

Press the black toggle switch on the right side of the Sea Salt Resurfacing Control Panel to the 'ON' position (I). The green light will illuminate when the Salt Flow is 'ON'.



TREATMENT AREA	INTENSITY LEVEL*	NUMBER OF PASSES
FACE, NECK & CHEST	3-5	2-3
PERI-ORBITAL	3	1
ARMS & HANDS	4-6	2-3
LEGS & TOPS OF FEET	4-6	3
BACK & SHOULDERS	4-6	3-4
ABDOMEN	4-6	3-4



STEP 2: AESTHETIC ULTRASOUND

1. **APPLY** 3-4 pumps of the SaltFacial DermMasque or preferred topical.
2. **APPLY** a thin layer of SaltFacial Gliding Gel (do not rub in).
3. **REMOVE** Aesthetic Ultrasound from its holder and remove blue cap.
4. **TURN** Aesthetic Ultrasound Module ON by pressing POWER button until the therapy screen turns on.

DO NOT USE FINGERNAILS TO PRESS BUTTONS AS THEY CAN PERMANENTLY DAMAGE THE MEMBRANE

5. **CONFIRM** treatment settings: 50% Duty Cycle / 8:00 minute Treatment Time / 1.85 W/cm²
6. **PLACE** ultrasound head on treatment area and move in slow circular motions until timer reaches 0:00 and treatment is complete.


To pause treatment, remove ultrasound head from treatment area. Treatment timer will also pause until contact has been re-established.

If contact has not been re-established within 3 minutes the unit will automatically power off (resetting the treatment timer).

7. **REMOVE** excess gel from treatment area with dry soft gauze or a soft cloth.

ROUTINE CLEANING: Clean transducer head after use and before putting protective cap over transducer.

NOTE: For sterile usage or when blood touches transducer, a sterile sleeve is recommended.



VISUAL & AUDIBLE CONFIRMATION OF CONTACT

When in full contact the blue light at the top of the Aesthetic Ultrasound Control Panel will illuminate, the treatment timer will automatically count down from 8:00, and there will be an audible beep every few seconds.

If contact is lost during treatment, apply more Gliding Gel and reposition the ultrasound head on the treatment area.



STEP 3: LED PHOTOTHERAPY

1. PLACE PROTECTIVE GOGGLES ON PATIENT

2. POSITION LED PANELS

Position LED Panels 4-6" from treatment area.

3. TURN LED PHOTOTHERAPY MODULE ON

Press Power button on far left side of control panel to turn power on.

4. SELECT TREATMENT MODE

Press Mode to select desired treatment mode.

Treatment mode will illuminate when selected and timer will automatically set


5. START TREATMENT

Press and hold On/Pause until white light on far right of LED Phototherapy Control Panel illuminates then immediately release.

To pause treatment press On/Pause. Press On/Pause again to resume treatment.



LED PHOTOTHERAPY: WAVELENGTH CHART



LED PHOTOTHERAPY: TREATMENT MODES

ACNE LIGHT THERAPY

WAVELENGTHS: BLUE + RED

TREATMENT TIME: 20 MINUTES

USE: As final step of SaltFacial Skin Renewal Therapy, or immediately post procedure (laser, microneedling, chemical peel, etc) to reduce redness and inflammation, stimulate cellular turnover, increase circulation, and promote healing.

SKIN TYPES: All***

PHOTODYNAMIC THERAPY

WAVELENGTHS: BLUE

TREATMENT TIME: 9 MINUTES

USE: As the final step of SaltFacial Skin Renewal Therapy, immediately prior to or post procedure (laser, microneedling, chemical peel, etc.) to destroy surface and sub-surface bacteria.

SKIN TYPES: All

COLLAGEN RESTORATION

WAVELENGTHS: RED + NEAR INFRARED

TREATMENT TIME: 20 MINUTES

USE: As final step of SaltFacial Skin Renewal Therapy, or immediately post procedure (laser, microneedling, chemical peel, etc) to reduce redness and inflammation, stimulate cellular turnover, increase circulation, and promote healing.

SKIN TYPES: All*

***CONTAINDICATED
FOR USE ON PATIENTS
WITH ACTIVE ROSACEA**

SKIN REJUVENATION

WAVELENGTHS: RED + AMBER


TREATMENT TIME: 20 MINUTES

USE: As final step of SaltFacial Skin Renewal Therapy, or immediately post procedure (laser, microneedling, chemical peel, etc) to reduce redness and inflammation, stimulate cellular turnover, and promote healing.

SKIN TYPES: All



CARE OF RECHARGEABLE DEHUMIDIFIER




Store extra fully-charged dehumidifier in a plastic bag and in a dry place.
NEVER LEAVE DEHUMIDIFIER PLUGGED IN FOR MORE THAN 15 HOURS



CARE OF THE LED ARTICULATING ARM & PANELS

- It is important to move the articulating arm carefully and NOT with force.
- Always support articulating arm between **JOINT 1** & **2** when loosening and tightening the LED assembly.
- Always loosen **JOINT 2** (middle elbow) before lowering LED arm.
- May cause damage to piston if not loosened - may void warranty.
- Do not rotate LED cable around **JOINT 1** or **3**.



CAUTION DROP HAZARD:

When loosened, LED head can cause injury or damage to system.



DAILY MAINTENANCE

1. **CLEAN** all accessory items used (hand piece, ultrasound head, LED goggles)

2. PURGE SYSTEM

OPEN the Salt Chamber door on the front panel and remove the salt canister from the holder - turn over the canister so the black Salt Cap faces up. Hold the canister lower than the level of tubing visible within the chamber.

TURN ON both the Salt Flow and System and set to Level 10.

HOLD your thumb over the opening on the hand piece until the system begins to 'click' and/or all visible signs of salt are cleared from the line.

TURN the Salt Flow and System OFF.


REPLACE the Salt Canister to the cradle by turning it over and gently seating it back in place. Some Salt will fall back into the tubing.

CLOSE the Salt Chamber door and the process is complete.

* Take this opportunity to check dehumidifier to make sure it is orange, not green

3. CHECK Rechargeable Dehumidifier

Confirm beads are **ORANGE**. If beads are **GREEN**, remove from system and plug into wall to recharge (See 'Care of Rechargeable Dehumidifier' for additional instructions).



WEEKLY MAINTENANCE

1. CLEAN NO-VAC EDUCTOR

REMOVE the No-Vac Eductor from The SaltFacial® System.

DIP small brush provided with The SaltFacial® System in alcohol and clean the 3 openings.

Note: Only use provided brush to clean No-Vac Eductor or hand piece. Do not use bottle brush or other brushes. Bottle brush or other brushes could damage equipment.

RUN under hot water for 1 minute, through all 3 openings. Use the provided small brush to clean all 3 openings.

Let dry overnight. Confirm no moisture remains in the No-Vac Eductor before using.

BEFORE reusing, confirm no moisture remains: Turn on the Salt Exfoliation toggle switch and thoroughly dry the No-Vac Eductor with the dry air from the No-Vac Air Connector Tube before reconnecting.

2. REPLACE RECHARGEABLE DEHUMIDIFIER

Replace rechargeable dehumidifier only if needed (ie. beads are GREEN).



CONTACT INFORMATION

TRAINING & SUPPORT CENTER

- Visit www.TheSaltFacial.com and click **LOGIN**
- Video Training Library
- Set Up/ Install Videos
- Procedure Videos
- Maintenance Videos

TECHNICAL SERVICE

- Visit www.thesaltfacial.com and click **LOGIN**
- Access the Troubleshooting Decision Tree

FOR REORDERING SUPPLIES OR MARKETING SUPPORT

(760) 334.5267 • info@saltfacial.com

